

Experiment Number-02

Name of the Experiment: To calibrate a Polarimeter and hence to determine the specific rotation of a sugar solution by means of a Polarimeter.

Data Collection:

Table for angular rotation

| Strength of sugar solution (%) | No. of obs | First reading with water (P) | Second reading with solution (Q) | Angular rotation (Q~P) | Mean angular rotation | Specific rotation (s) |
|---------------------------------------|-------------------|-------------------------------------|---|-------------------------------|------------------------------|------------------------------|
| 20% | | | | | | |
| 10% | | | | | | |
| 5% | | | | | | |

$$\text{Specific rotation} = \frac{10\theta}{lc}$$

Experiment Number-3

Experiment Name: Determination of moment of inertia of a flywheel about its axis of rotation.

Data Collection:

Table 1: Determination of n_1 , n_2 and t :

| Mass M gm | Height h cm | No. of Revolutions n_1 | Average n_1 | No. of Revolution s n_2 | Average n_2 | Time t | Average t | Moment of Inertia I | Average I |
|-----------------|-------------------|--------------------------------|------------------|---------------------------------|------------------|-----------|--------------|---------------------------|--------------|
| | | | | | | | | | |
| | | | | | | | | | |
| | | | | | | | | | |
| | | | | | | | | | |
| | | | | | | | | | |
| | | | | | | | | | |
| | | | | | | | | | |
| | | | | | | | | | |
| | | | | | | | | | |

Table 2: Determination of the radius of the axle:

| LSR (X) (cm) | VSD | VC (cm) | VSR Y =(VSD×VC) (cm) | Total Reading = X+Y cm | Average Diameter cm | Radius r cm |
|--------------------|-----|------------|----------------------------|------------------------------|---------------------------|-------------------|
| | | | | | | |
| | | | | | | |
| | | | | | | |

Calculation:

$$\omega = \frac{4\pi n_2}{t} \text{ and } I = \frac{2mgh - m \omega^2 r^2}{\omega^2 (1 + \frac{n_1}{n_2})}$$

Experiment Number-4

Name of the Experiment: To determine the Specific resistance of a wire using a Meter Bridge

Data Collection:

Table 1: Reading for the Galvanometer balance point.

| No. of observation | Value of resistance, R(Ω) | Length, l cm | (100 - l) cm | $X = \frac{R(100-l)}{l} (\Omega)$ | Mean X(Ω) |
|--------------------|------------------------------------|--------------|--------------|-----------------------------------|--------------------|
| 1. | | | | | |
| 2. | | | | | |
| 3. | | | | | |
| 4. | | | | | |
| 5. | | | | | |
| 6. | | | | | |
| 7. | | | | | |
| 8. | | | | | |
| 9. | | | | | |
| 10. | | | | | |

Table 2: Reading for the radius of the experimental wire. (Using Slide Calipers)

| No. of obs. | Liner scale reading (L.S.R) cm | Circular scale divisions (C.S.D) | Least count (L.C) cm | Circular scale reading (C.S.R) = (CSD x L.C) cm | Total diameter D cm | Mean diameter D cm | Mean radius r=D/2 cm |
|-------------|--------------------------------|----------------------------------|----------------------|---|---------------------|--------------------|----------------------|
| 1 | | | | | | | |
| 2 | | | | | | | |
| 3 | | | | | | | |

Specific Resistance

$$\rho = \frac{\pi r^2 X}{L}$$

Experiment Number-5

Name of the Experiment: To determine the e.m.f of a cell with a potentiometer of known resistance

Data collection:

- (A) Resistance of the potentiometer wire, $R = 72 \Omega$
(B) Total length of the potentiometer, $L = 1000 \text{ cm}$.

| No. of obs. | Miliammeter readings = i mA | Null points | | Total length for balance l cm | e.m.f of the cell $E = \frac{iRl}{1000L}$ Volts. | Mean E' Volts. |
|-------------|-----------------------------------|----------------|---------------------|---------------------------------------|--|------------------------|
| | | On wire number | Scale reading in cm | | | |
| 1 | | 10th | | | | |
| 2 | | 9th | | | | |
| 3 | | 8th | | | | |
| 4 | | 7th | | | | |

Experiment Number-6

Name of The Experiment: To compare the e.m.f of two cells with a potentiometer

Data collection:

Table 1: Data for the value of E_1/E_2

| No. of observations | Cell | Null points | | Total length in cm | $E_1/E_2 = l_1/l_2$ | Mean E_1/E_2 |
|---------------------|------------------|----------------|---------------------|--------------------|---------------------|----------------|
| | | On wire number | Scale reading in cm | | | |
| 1 | First (E_1) | | | | | |
| | Second (E_2) | | | | | |
| 2 | First (E_1) | | | | | |
| | Second (E_2) | | | | | |
| 3 | First (E_1) | | | | | |
| | Second (E_2) | | | | | |

Experiment Number-7

Name of The Experiment: To determine the wavelength of various spectral lines by a spectrometer using a plane diffraction grating.

Data Collection:

Table 1: Data for determination of grating constant

[illegible]

Table 2: Data for determination of wavelength of spectral lines for Mercury discharge tube

[illegible]

Experiment Number-09

Name of the Experiment: To determine the value of unknown resistance and verify the laws of series and parallel resistance by Post Office Box.

Data collection:

Table 1: Reading for the Unknown resistance.

| Value of P (Ω) | Value of Q (Ω) | Value of R (Ω) | Value of S $= \frac{Q \times R}{P}(\Omega)$ | Mean, S(Ω) |
|-------------------------|-------------------------|-------------------------|--|------------------------|
| 10 | 10 | | | |
| 10 | 100 | | | |
| 100 | 1000 | | | |

Table 2: Reading for the series resistance, S_s

| Value of P (Ω) | Value of Q (Ω) | Value of R (Ω) | Series resistance, $S_s = \frac{Q \times R}{P}(\Omega)$ | Mean $S_s(\Omega)$ |
|-------------------------|-------------------------|-------------------------|--|-----------------------|
| 10 | 10 | | | |
| 10 | 100 | | | |
| 100 | 1000 | | | |

Table 3: Reading for the parallel resistance, S

| Value of P (Ω) | Value of Q (Ω) | Value of R (Ω) | Parallel resistance $S_p = \frac{Q \times R}{P}(\Omega)$ | Mean, $S_p(\Omega)$ |
|-------------------------|-------------------------|-------------------------|---|------------------------|
| 10 | 10 | | | |
| 10 | 100 | | | |
| 100 | 1000 | | | |

